



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Non-pressure tank for storing and transporting liquid nitrogen AIR LIQUIDE TR26 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in France,
- **Zbiornik TR26** is a non-pressurized vessel designed for storing and transporting liquid nitrogen. Their lightweight aluminum construction with a composite plastic bottleneck ensures the vessel's light weight, resistance and minimal nitrogen evaporation. Polyurethane paint coating ensures high quality and prolongs the life of the vessel.
- **Specification:**
  - Made of aluminium,
  - Productive capacity: 26 l,
  - Bottleneck diameter: 50 mm,
  - Daily evaporation: 0.2 l,
  - Static holding time: 130 days,
- Has visible marks of usage,
- Mobile on a mobile base with 5 wheels including 2 with a brake,
- **The set includes:**
  - Mobile base,
- Dimensions:
  - Diameter: 39 cm,
  - Height: 67 cm,
- Weight: 20 kg,
- Cat 02,
- Provided with a valid Technical Passport,
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

