



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Non-pressure tank for storing and transporting liquid nitrogen AIR LIQUIDE TR26 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in France,
- **Zbiornik TR26** is a non-pressurized vessel designed for storing and transporting liquid nitrogen. Their lightweight aluminum construction with a composite plastic bottleneck ensures the vessel's light weight, resistance and minimal nitrogen evaporation. Polyurethane paint coating ensures high quality and prolongs the life of the vessel.
- **Specification:**
 - Made of aluminium,
 - Productive capacity: 26 l,
 - Bottleneck diameter: 50 mm,
 - Daily evaporation: 0.2 l,
 - Static holding time: 130 days,
- Has visible marks of usage,
- Mobile on a mobile base with 5 wheels including 2 with a brake,
- **The set includes:**
 - Mobile base,
- Dimensions:
 - Diameter: 39 cm,
 - Height: 67 cm,
- Weight: 20 kg,
- Cat 02,
- Provided with a valid Technical Passport,
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

