



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



A tank for long-term biological samples storage AIR LIQUIDE GT35 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in France,
- **GT35 tank** for long-term storage of biological samples in straws or vials is characterized by very low liquid nitrogen consumption. Dedicated to storing a large number of samples in a small volume,
- **Specification:**
 - Made of aluminium,
 - Productive capacity: 33,6 l,
 - Bottleneck diameter: 50 mm,
 - Daily evaporation: 0,09 l,
 - Static holding time: 350 days,
 - Containers' number: 6,
 - Containers' diameter: 38 mm,
 - Containers' height: 280 mm,
 - The number of levels: 1,
 - Maximum capacity of 0.25 ml straws: 3120,
 - Maximum capacity of 0.5 ml straws: 3120,
 - Max. capacity of 2 ml vials: 180,
- Has visible marks of usage,
- Mobile on a mobile base with 5 wheels including 2 with a brake,
- **The set includes:**
 - 6 x 280 mm containers,
 - Mobile base,
- Dimensions:
 - Diameter 54 cm,
 - Height: 90 cm,
- Weight:
 - An empty container: 15 kg,
 - A full container: 43 kg,
 - Mobile trolley: 6 kg,
- Cat 02,
- Provided with a valid Technical Passport,
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,



