



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MELAG VACUKLAV 31 autoclave (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany (2000),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2500 W,
- **MELAG Vacuklav 31** is an efficient steam autoclave offering 4 work programs and 2 test programs,
- **Product's characteristics:**
  - B class,
  - Automatic,
  - Pre-vacuum and fractional vacuum generated by a vacuum pump,
  - Additional vacuum drying is possible,
  - Preheat option helps maintain temperature between cycles,
  - Equipped with indicator lights to indicate unit errors,
  - 2.5" LCD display,
  - Foiled control panel,
  - On the lid of the autoclave, a place to put down the tray,
  - Water drain connector on the front of the unit,
  - Sockets:
    - For Melaprint printer **(not included)**,
    - For Meladem ion exchange system **(not included)**,
    - 230 V outlet,
  - Made of stainless steel materials,
  - Stable on legs, including front legs with height regulation,
  - German keyboard on the control panel,
- **It has 4 work programs and 2 test programs:**
  - Unpacked tools: 2 bar pressure - time 4 min - temperature 134 °C,
  - Packed tools: 2 bar pressure - time 5 min - temperature 134 °C,
  - Textiles: 1 bar pressure - time 17 min - temperature 121 °C,
  - Disinfection: 0.3 bar pressure - time 8 min - temperature 105 °C,
  - The Bowie & Dick test: vapor permeation test to check the efficiency of the autoclave by the user,
  - Leakage test: vacuum test to check the leakage of the autoclave,
- **Specification:**
  - Sterilization chamber capacity: 14,5 L,
  - Maximum load capacity: 4 kg of instruments including trays or 1.5 kg of textiles,
  - Water consumption during sterilization: 300 - 400 ml, depending on the program and loading level,
- **In the set:**
  - English instruction manual (PDF),
  - Trays of dimensions: 18 x 28 x 2 cm - 2 pieces,
  - Trays' handles,
  - Stand with a handle for 3 trays,
  - Brand new condensation tank,
- Dimensions:
  - Of the autoclave: 63 x 45 x 45 cm,
  - Of the chamber: Ø 23 cm x 34 cm,
- Weight: 44 kg,
- Cat 03,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport and ready to work,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

