



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



JORG & SOHN EXAFLEX 3000 operating table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- **JORG & SOHN EXAFLEX 3000** is a mobile multifunctional operating table designed for general use in hospitals and ambulatory care units,
- **Medically certified (memory visco) upholstery in black (can be changed to new one for extra fee),**
- 5 - segment,
- **Manual regulation:**
 - Height with the use of hydraulic pump lever in the range of: 78 - 108 cm,
 - Trendelenburg and anti-Trendelenburg positions by using a gas actuator in the range of: $\pm 33^\circ$,
 - Lateral tilt using crank mechanism in the range of: 18° ,
 - Backrest by gas actuator in the range: from -30° to 80° ,
 - Headrest in multiple surfaces,
- **It has:**
 - Ability to mount x-ray cassettes,
 - Rails for additional accessories mounting,
 - Central brake,
- The base is made of chrome-nickel steel,
- **It is possible to equip the table (for an extra charge) with additional accessories, examples shown in the last photos,**
- Dimensions:
 - Of the bedspace: 200 x 54 cm,
 - Of the base: 103 x 49 cm,
- Maximum load: 160 kg,
- Weight: 190 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

