



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **JORG & SOHN EXAFLEX 3000 operating table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- **JORG & SOHN EXAFLEX 3000** is a mobile multifunctional operating table designed for general use in hospitals and ambulatory care units,
- **Medically certified (memory visco) upholstery in black (can be changed to new one for extra fee),**
- 5 – segment,
- **Manual regulation:**
  - Height with the use of hydraulic pump lever in the range of: 78 - 108 cm,
  - Trendelenburg and anti-Trendelenburg positions by using a gas actuator in the range of:  $\pm 33^\circ$ ,
  - Lateral tilt using crank mechanism in the range of:  $18^\circ$ ,
  - Backrest by gas actuator in the range: from  $-30^\circ$  to  $80^\circ$ ,
  - Headrest in multiple surfaces,
- **It has:**
  - Ability to mount x-ray cassettes,
  - Rails for additional accessories mounting,
  - Central brake,
- The base is made of chrome-nickel steel,
- **It is possible to equip the table (for an extra charge) with additional accessories, examples shown in the last photos,**
- Dimensions:
  - Of the bedspace: 200 x 54 cm,
  - Of the base: 103 x 49 cm,
- Maximum load: 160 kg,
- Weight: 190 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

