



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ALVO SONATA 4-03 operating table. (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Poland (2012),
- Power supply: 230 V and via 2 in-built 12V 7.2 Ah 20HR batteries,
- Frequency: 50/60 Hz,
- **ALVO Sonata** table systems was developed to optimize workflow in the operating room. The system allows to swap tops - removing and attaching the top allows to use a cart for transporting the top, which is compatible with a wide range of bases (flat, mobile, rotating). SONATA system tables are designed to support the patient during all surgical procedures and operations: in general surgery, radiology, endoscopy, neurosurgery, gynecology, urology and orthopedics. The regulation of the table's functional movements is carried out by a reliable electro-hydraulic system. Functional tabletop movements are supported by gas springs or mechanically,
- **Medically certified upholstery in black with an antibacterial BioCote layer,**
- 5 - segment,
- Electrohydraulic drive,
- **Electrically regulated with the use of a wired remote control in many surfaces:**
 - Height in the range of: 112 - 80 cm,
 - Trendelenburg and anti-Trendelenburg positions in range: 40°,
 - Lateral tilt in range: 30°,
 - "0" position button,
- **Manual regulation of:**
 - Backrest angle in the range: from -45° to 85°,
 - Footrest angle in the range: from -90° to 20°,
 - Angle of footrests gape in the range: 180°,
 - Angle of the headrest in the range of: -45° to 45°,
 - Renal bench elevation: 15 cm,
 - Longitudinal movement: 30 cm,
- **It has:**
 - X-RAY transparency,
 - Rails for additional accessories attachment,
 - Ability to mount X-RAY cassettes,
 - An in-built charger,
- Operation type: intermittent 2/18 min,
- Made of stainless, acid-resistant steel,
- Mobile with a central brake,
- **Possibility to equip the table (for an extra charge) with additional accessories, examples shown in the last photos,**
- **In the set:**
 - 2 x transportation cart for **ALVO 4-09-006** top,
 - Second replacement **ALVO 4-03-5** top,
 - Wired remote control,
 - Power cord,
- Dimensions:
 - Of the bedspace: 210 x 50 cm,
 - Of the base: 57 x 99 cm,
- Maximum load: 275 kg,
- Weight: 250 kg,
- Cat 02,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport and ready to work,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

