



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **J&J CIDEX 2027 sterilization container (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- **J&J Cidex** is a container for sterilization and disinfection of medical instruments,
- Can be used for manual cleaning or high level disinfection,
- Temperature resistant up to 121°C,
- Reduces the evaporation of steam and vapors due to the unique rounded edge,
- Easy cleaning thanks to the wide radius corner
- Capacity: 14.1 L,
- Blue,
- **In the set:**
  - The cover,
- Dimensions:
  - Internal: 14 x 66,5 x 13,4 cm,
  - Overall: 19 x 72,3 x 13,8
- Weight: 1,5 kg,
- Cat 01,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport and ready to work,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

