



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **IVAC P1000 infusion pump - syringe (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- **Real photos of the product,**
- Made in England,
- Power supply: 230 V and via 1 new, in-built 6 V / 2,8 Ah battery,
- Frequency: 50/60 Hz,
- Power: 12 VA,
- **IVAC P1000** is a syringe pump designed for precise administration of continuous infusions. Extensive alarms and sophisticated monitoring systems, combined with simple operation, make the device suitable for general and critical infusions,
- **The pump has:**
  - Foiled control panel and liquid crystal display,
  - German menu,
  - Hanging handle,
  - Comfortable transportation handle,
  - In-built audible alarms,
- **Adjusted for following syringes:**
  - **10 ml, 20 ml, 30 ml, 50 ml:**
    - DB PLASTIPAK,
    - TERUMO,
    - B BRAUN OMNIFIX,
    - B BRAUN PERFUSOR,
    - SHERWOOD MONOJECT,
    - RR PRONTO,
  - **50 ml:**
    - FRESENIUS INJECTOMAT,
    - ONCE,
- **In the set:**
  - German instruction manual (PDF),
  - Power cord,
- Overall dimensions: 24 x 37 x 10,5 cm,
- Weight: 3 kg,
- Cat 02,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport and ready to work,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

