



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED L-120 negatoscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 65 W,
- **Specification:**
 - Single - frame,
 - Equipped with 5 fluorescent lamps: (L 13W/25),
 - Four-stage light intensity control,
 - Features horizontal and vertical blinds,
 - Maximum frame size: 41 x 42 cm,
 - Negatoscope for mounting on the wall,
- Overall dimensions: 14 x 61 x 66 cm,
- Weight: 18,5 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

