



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MERIVAARA OP1650 operating table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Finland,
- **MERIVAARA OP1650 operating table** is an easy-to-use, manual operating table designed for many surgical disciplines,
- **New upholstery (memory visco) with medical certificate in black,**
- 4 - segment,
- **Manually regulated in many surfaces:**
 - Height in the range: 73 - 100 cm,
 - Headrest angle in the range: from -45° to +25°,
 - Backrest angle in the range: from -15° to +90°,
 - Lateral tilt: ±20°,
 - Trendelenburg position: 25°,
 - Anti-Trendelenburg position: 18°,
 - Leg section in the range: from -90° to 0°,
- It has:
 - X-Ray transparency,
 - Rails for additional equipment attachment,
- The base made of stainless steel,
- Mobile with a brake,
- **The table can be equipped with extra accessories (at additional cost); examples visible in the last photos,**
- Dimensions:
 - Of the bedspace: 195 x 50 cm,
 - Of the base: 62 x 106 cm,
- Maximum load: 180 kg,
- Weight: 135 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

