



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SECA CT8000P ECG device (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Italy,
- Power supply: 230 V,
- Emergency power supply: via 1 new, in-built 12V / 2,3 Ah battery,
- Frequency: 50/60 Hz,
- **SECA CT8000P** is a high-quality interpretive ECG with LCD screen and alphanumeric keyboard. Ideal for use in clinics and hospitals with large numbers of patients who can rely on a fully functional ECG.
- **Specification:**
 - 12 - lead,
 - Features paperless ECG options in connection with seca Passport ECG software - compatible with leading EMR, emis, SystemOne, vision, isoft and other systems,
 - One-touch automatic or instant real-time operation,
 - Easy trace selection from 12 x 1, 6 x 2 or popular 4 x 3 plus Rythm Strip one-page report,
 - Features high quality digital filters for clear traces (grid, baseline and muscle tremor),
 - Internal memory (max 45 records),
 - All parameters, including filters, can be easily programmed via the on-screen menu for one-touch ECG,
 - SVGA output for connecting an external color monitor,
 - Dedicated one-touch keypads for storage, printing and transfer,
 - Easy-to-load paper tray (220 A4 pages),
- **It has normal signs of use including 2 very small flaws on the screen that do not have a major impact on the comfort of operating the device and reading the results,**
- **Menu to be chosen in: EN, DE, FR, US, IT, ES, PR, S, NL, HU, DK, FN, SE,**
- **English keyboard on the control panel,**
- **In the set:**
 - Instruction manual in EN, DE, FR, ES, PO (PDF),
 - 10 - lead ECG cable,
 - 10 x jaw clamp,
 - 90 x disposable electrodes,
 - Power cord,
 - Original transport bag,
- Dimensions: 33 x 38 x 8 cm,
- Weight: 6 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

