



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## LED CURING LIGHT UNIT DB686 MOCHA polymerization lamp (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Power supply via new, in-built DC 3,7 V 2200 mAh battery,
- Frequency: 50/60 Hz,
- **DB-686 MOCHA** - is a wireless lamp with a modern design,
  - Wavelength: 420 nm - 480 nm,
  - Light power: 1200 mW/cm<sup>2</sup>,
  - LED power: 8 W max,
  - High power LED, pure dental blue light,
  - Digital LCD display,
  - Built-in soft start,
  - Regulated run time,
  - Low battery detection,
  - Long standby time,
  - Overheat protection,
  - Light source detection,
- **In the set:**
  - Power adapter,
- Dimensions
  - Overall: 12,5 x 11,5 x 27,5 cm,
  - Of the carton box: 19 x 24,3 x 9 cm,
- Set's weight: 2 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

