



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ELEMENT 6 diagnostic-treatment table (New)

- A brand new product,
- Made in Poland,
- ELEMENT 06 stationary rehabilitation tables feature an ABS plastic case, which gives them a modern look and ensures easy maintenance,
- Standard features include a four-part recliner with a wide range of regulation by gas springs, made of plywood,
- Height regulation of the table is performed by an electric actuator,
- The table has a stable base and fully ball-bearing regulation mechanisms,
- There is a free option to order upholstery in the color of your choice, except for "PRESTIGE" colors which are available at extra cost,
- As additional options there are available (at extra cost):
 - Manually regulated armrests,
 - Lowered side shelves,
 - Hanging pad,
 - Mobile wheels system,
 - Foot remote control,
 - Frame control,
- Technical parameters:
 - Length 200 cm,
 - Width 65/70/80 cm,
 - Height: 55 - 106 cm,
 - Headrest angle regulation: -75° to 45°,
 - Breast section angle regulation: -18° to 80°,
 - Lumbar section break angle regulation: -20° to 84°,
 - Footrest angle regulation: -45° to 50°,
 - Trendelenburg position: 0° to 20°,
 - Anti-Trendelenburg position: -18° to 0°,
 - Own weight 113 kg,
 - Maximum load 160 kg,
- It has a valid Technical Passport issued,
- 24 months' guarantee for the construction,
- 12 months' guarantee for hydraulic accessories,

