



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ELEMENT CH diagnostic-treatment table (New)

- A brand new product,
- Made in Poland,
- ELEMENT 06 stationary rehabilitation tables feature an ABS plastic case, which gives them a modern look and ensures easy maintenance,
- Standard features include a six-section recliner, made of plywood,
- Height regulation of the table is performed by an electric actuator,
- The advantage of the table is the possibility to lift with gas springs, separately or together the lumbar and chest sections,
- An important feature in the offered model is the wide range of regulation of the headrest,
- The modern solution of upward and downward position regulation with tilt angle adjustment, allows comfortable positioning of the patient,
- The table has a stable base and fully ball-bearing regulation mechanisms,
- There is a free option to order upholstery in the color of your choice, except for "PRESTIGE" colors which are available at extra cost,
- As additional options there are available (at extra cost):
 - Manually regulated armrests,
 - Pad hanger,
 - Mobile wheels system,
 - Foot remote control,
 - Frame control,
- Technical parameters:
 - Length 195 cm,
 - Width 65/70 cm,
 - Height: 48 - 100 cm,
 - Headrest angle regulation: -30° to 30°,
 - Headrest height regulation: 26 cm,
 - Headrest tilt: 6 to 12 cm,
 - Backrest angle regulation: 0° to 15°,
 - Mid-section angle regulation: 0° to 25°,
 - Own weight 90 kg,
 - Maximum load 150 kg,
- It has a valid Technical Passport issued,
- 24 months' guarantee for the construction,
- 12 months' guarantee for hydraulic accessories,

