



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **NGP 11 WS negatoscope (New)**

- A brand new product,
- Made in Poland,
- The NGP-11 WS wall-mounted negatoscope is a highly functional and damage-resistant device,
- It emits daylight, completely friendly to the radiologist's eyes,
- **Advantages of using the NGP 11 WS negatoscope:**
  - They can be used in both operating rooms and treatment rooms,
  - It is a completely easy to disinfect, sealed product,
  - It is made of acid-resistant sheet metal,
  - They have all the required certificates,
- **Technical data:**
  - Supply voltage: 230 V 50 Hz,
  - Power consumption: 65 VA,
  - Screen luminance: 6000 cd/m<sup>2</sup> ±15%,
  - Luminance regulation: YES, smooth 10-100%,
  - Screen dimensions: 36 x 43 cm,
  - Number of frames: 1,
  - Protection class: I,
  - Internal dimensions: 435 x 545 x 50 mm,
  - Dimensions: 470 x 580 x 74 mm,
  - Light source: LED panel,
  - Weight: 6 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

