



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NGP 31 WS negatoscope (New)

- A brand new product,
- Made in Poland,
- The NGP-31 WS wall-mounted negatoscope is a highly functional and damage-resistant device,
- It emits daylight, completely friendly to the radiologist's eyes,
- The uniform, evenly illuminated screen will enable accurate diagnosis of the patient's disease thanks to clear contrasts,
- **Advantages of using the NGP 31 WS negatoscope:**
 - They can be used in both, operating and treatment rooms,
 - It is a completely easy to disinfect, sealed product,
 - It is made of acid-resistant sheet metal,
 - They have all the required certificates,
- **Technical data:**
 - Supply voltage: 230 V 50 Hz,
 - Power consumption: 180 VA,
 - Screen luminance: 6000 cd/m² ±15%,
 - Luminance regulation: YES, smooth 10-100%,
 - Screen dimensions: 108 x 43 cm,
 - Number of frames: 3,
 - Protection class: I,
 - Internal dimensions: 1155x 545 x 50 mm,
 - Dimensions: 1190 x 580 x 74 mm,
 - Light source: LED panel
 - Weight: 14,0 kg,
- Provided with an up-to-date Technical Passport,
- Guarantee: 12 months,

