



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ENTERMED ERGOMED 600 laryngological unit (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the Netherlands,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2300 VA,
- The **ENTERMED ERGOMED 600 unit** uses top-quality materials, including a wooden segment case and a top made of luxurious MARLAN marble. As a result, the unit has an elegant design and is very efficient. The unit's construction is modular. It consists of two standing segments and two wall-mounted cabinet segments.
- **Segment I:**
  - Foiled control panel,
  - Two level top,
  - Two pull-out shelves on slides located under the upper top, configurable aluminum trays of customer's choice,
  - The drawer is double tier with dividers, ability to configure the dividers as desired, plus a sight glass rack,
  - 2 drawers with opening for used instruments and waste container,
  - Mirror warmer with automatic shut-off after 3 seconds,
  - Rigid optics heater,
  - Integrated halogen light source automatically switched on by raising the light guide
  - Small basin activated via control panel,
  - Slipper - suction system with 1.5l disposable cartridge jar (Baxter).
  - Has 2 power outlets for connecting additional accessories,
  - The size of the segment : 144 x 80cm the distance of the top from the floor 94 cm,
- **Segment II :**
  - Sanitary area with a washbasin,
  - Two level countertop,
  - The washbasin with a diameter of 38 cm and a depth of 14 cm with a faucet for water,
  - Pull-out shelf on slides located under the upper countertop, aluminum trays can be configured as desired,
  - Cabinet with hinged doors with a metal basket installed, dimensions : 24 x 33 cm and depth of 40 cm,
  - It has four drawers with dimensions:
    - 19 x 46 cm and 7 cm deep (3 pcs).
    - 19 x 46 cm and 18.5 cm deep,
  - Independent light source Entermmed HL-150 with smooth light intensity regulation,
  - Possibility to connect a headlamp automatically switched on by raising the lamp (no lamp included),
  - Segment size : 123 x 55.5 cm top to floor distance 94 cm,
- **Segment III hanging** (with an original hole for a monitor or additional hardware):
  - Possibility to change the door panel with a hole (where the monitor is now, to a door with hinges (for additional fee),
  - The monitor can be mounted on a fixed mast or on a fixed mast with a movable arm coming out of the tabletop (for additional fee)
  - A shelf for an additional light source is located on the underside with an attachment,
  - Illumination under the cabinet (two lamps with OSRAM Linestra 230V 35W bulbs)
  - The segment has two cabinets with shelves (the height of the shelves can be regulated)
  - Dimensions of the cabinet compartment : 33 x 56 cm and depth of 32 cm,
  - Segment's dimensions : 130,5 x 57 x 67,5 cm,
- **Segment IV hanging with a negatoscope:**
  - Negatoscope by ENTERMED:
  - Made in the Netherlands,
  - 230V power supply,
  - Equipped with six fluorescent lamps of 15 W each,
  - It has a cabinet with a shelf with a chamber size of 28 x 56 cm and depth of 32 cm (the height of the shelf can be regulated),
  - Illumination under the cabinet (two lamps with OSRAM Linestra 230V 35W bulbs),
  - Segment's dimensions: 123 x 37,5 x 67,5 cm,
- **The unit can be extra equipped with the following (for additional fee):**
  - Video camera processor,
  - Video camera head,
  - Endoscope,
  - Medical monitor,
  - Head lamp,
  - Operating microscope,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

