



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHMITZ MEDI-MATIC operating table 125.200 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German manufacture,
Upholstery with a medical certificate in black (upholstery can be replaced with a new one, at an additional cost),
5 - segmented,
Operating table
SCHMITZ MEDI-MATIC 125.200
is an easy-to-use, manual operating table designed for many surgical fields,
Adjustable manually in many planes with hydraulic and gas cylinders:
Height in range: 70 - 100 cm,
Back rest angle in the range: from -15° to +60°,
Lateral tilt: $\pm 20^\circ$,
Trendelenburg position: 30°,
Anti-Trendelenburg position: 25°,
Headrest adjustment angle in range: 40°,
Leg segment in range: from -90° to +90°,
Lateral recoil in range: 90°,
Features:
X-Ray translucency,
Split leg segment laterally inverted,
Rails for attaching additional accessories,
Stainless steel base,
Mobile with central brake,
Possibility to retrofit the table (for an extra charge) with additional accessories, examples shown in the last photos,
Dimensions:
Bed: 200 x 50 cm,
General: 203 x 56 x 100 cm,
Weight: 180 kg,
Has up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

