



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MIRAE MRK-2000 Auto kerato - refractometer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Korea,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 53 W,
- Automatic refractometer with keratometry, has the necessary functions to perform full diagnostics,
- Technical specifications:
 - Refractometry:
 - Sphere: -25.0 dpt to +22.0 dpt (for VD=12mm, grading 0.12-0.25 dpt),
 - Cylinder: 0.00 - ±10.0 dpt (grading 0.12-0.25 dpt),
 - Keratometry:
 - Corneal power: 33.00 ~ 67.50 dpt (grading 0.05/0.12/0.25 dpt),
 - Memory for 10 measurement results for L/R eye,
 - Built-in thermal printer with 6 cm paper diameter,
- English language menu,
- Included:
 - Power cord,
- Dimensions: 27 x 50 x 47 cm,
- Weight: 20 kg,
- Has an up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

