



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **ASPEL B612 Exercise test treadmill (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2300 VA,
- ASPEL B612 treadmill is used for exercise testing and cardiac rehabilitation equipped with world unique solutions to improve medical practice,
- Technical specifications:
  - Cooperation with CardioTEST exercise testing system and AsTER cardiac rehabilitation system,
  - Belt speed adjustable from 0 to 20 km/h,
  - Tilt angle adjustable from 0% to 25%,
  - Usable belt length: 130 cm, width 500 cm,
  - Interlocked safety switch, stops the tape running,
  - ZERO START function,
  - Unique incline control system - built-in inclinometer to ensure automatic leveling of the treadmill and leveling of uneven ground,
  - Built-in ECG module with the following parameters:
    - 10 standard leads in diagnostic mode,
    - 2 leads in training mode,
    - Sensitivity 2.5/5/10/20 mm/mV,
    - Offset 25/50/100 mm/s,
    - Digital filters 25Hz, 35Hz, 50Hz,
    - Digital isoelectric line filter,
  - Control using computer via RS-232 port,
- Kit equipped with CARDIOTEST 2.71.00 software (2006),
- Ability to select any language in the system,
- Included:
  - User manual in Polish (PDF),
  - Power cord,
  - RS-232 cable,
  - 10-lead ECG cable,
  - Computer + operator monitor,
  - Operator's table,
- Dimensions:
  - Treadmill: 82 x 200 x 125 cm,
  - Operator station: 90 x 60 x 80 cm,
- Maximum load capacity: 180 kg,
- Weight: 170 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU

