



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **FISHER & PAYKEL MR850 Ventilator humidifier heater (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in New Zealand,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 150 W,
- One button selects the optimal temperature and humidity level for adult, pediatric and neonatal patients,
- Easy setup with simple controls and easy-to-read display,
- Alarm system,
- Advanced algorithms and flow detection capability,
- Technical specifications:
  - Invasive mode:
    - Chamber outlet: 35.5 - 42 °C / airway: 35 - 40 °C,
  - Non-invasive mode:
    - Chamber outlet : 31 - 36 °C / airway : 28 - 34 °C,
  - Temperature range:
    - 10 - 70 °C / accuracy of 0.3 °C (range 25 - 45 °C),
  - Lacks temperature sensor and breathing circuit heater adapter (can be retrofitted at additional cost),
- Included:
  - User manual in English (PDF),
  - Power cord,
- Dimensions:
  - Overall: 14 x 17 x 16 cm,
- Weight: 1.5 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

