



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **WELCH ALLYN PROPAQ CS 242 Patient monitor / Cardiac monitor (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in USA,
- Power supply: 230 V,
- Emergency power supply: through 2 new built-in 8 V / 2.7 Ah batteries,
- Frequency: 50/60 Hz,
- PROPAQ CS 242 is a compact patient monitor that allows the operator to control many of the patient's vital signs even under extreme conditions,
- The device has:
  - 8.5" LCD touchscreen display,
  - The ability to connect to a defibrillator via the DEFIB SYNCHRO jack,
  - Audible and visual alarms,
  - Handle for moving and hanging,
- Equipped with the following modules with full wiring and sensors:
  - ECG / ECG RESP - New 3 - lead cable with snap type connectors,
  - SPO2,
  - NIBP - New cuff with air line,
  - TEMP 1/2 - No sensor included (can be retrofitted at additional cost),
- English language menu,
- Included:
  - User manual in English (PDF),
  - Power cord,
  - 19 V power supply,
- Dimensions: 24 x 20 x 20 cm,
- Weight: 6 kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.



