



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



HETTICH MIKRO 220R Laboratory centrifuge (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- (YOM 2019),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 900 W,
- Foiled control panel with digital readable display,
- Automatic rotor recognition,
- High acceleration,
- Short spin cycles (pulse button),
- Rapid cooling,
- Brushless induction motor for years of maintenance-free operation,
- Technical specifications:
 - Infinitely variable speed control max. up to 18,000 rpm, with rotor used speed control from 500 to 6,000 rpm, max. RCF 4,146,
 - Max. centrifugal acceleration: 31514 x g,
 - Max. Capacity - angle rotor: 48x1.5/2.0 ml; 6x50 ml,
 - Range of generated temperature in the range: -20 / +40°C,
 - 9 acceleration and deceleration speeds,
 - Timed operation switch in range: 1 sec - 99 min, continuous run, short cycle (pulse key),
- Centrifuge equipped with 35° angle rotor HETTICH 1015:
 - Test tubes: 12 x 35 g,
 - 18 x 88 mm,
- Included:
- Instruction manual in English, German, French, Italian (PDF),
 - Power cord,
 - HETTICH 1015 rotor,
 - Inbus key,
- Dimensions: 33 x 65 x 31 cm,
- Weight: 42 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12 months for domestic and international markets for a surcharge,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

