



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DANTEC DYNAMICS DANFLOW 1000 Uroflowmeter (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Czech Republic,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Uroflometer connected by a 3m cable to a compact control and measurement module in a stainless steel housing,
- Designed to automatically measure the course of uroflowmetry,
- Prints the test on thermal paper (56mm),
- Measurement range is 0-50 ml/s, with accuracy of $\pm 2\%$,
- Standard stainless steel stand with urinal module,
- DanFlow 1000 can be used for male or female patients,
- Adjustable height,
- Included:
 - Power cord,
 - Weight,
 - Removable urinal,
 - Urine container with handle,
 - Control panel with built-in printer,
- Dimensions:
 - Uroflometer: 20 x 16 x 5 cm,
 - Urinal with instrumentation: 32 x 30 x 40 cm,
- Weight: 7 kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

