



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **CSM-2 DANMETER Cerebral state monitor (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the Netherlands,
- Power supply: 230 V,
- Emergency power supply: via 1 built-in battery (9 V),
- Frequency: 50/60 Hz,
- CSM Cerebral State Monitor records EEG using three EEG electrodes placed on the patient's forehead. The frequency content obtained from the EEG is converted into CSI (brain state indicator) within seconds using a neuro-fuzzy logic processor,
- Technical specifications:
  - EEG sensitivity:  $\pm 400 \mu\text{V}$ ,
  - CMRR:  $>140 \text{ dB}$ ,
  - CSI value: 100 linear Update every 1 s, filtering 6 - 42.5 Hz,
  - Sensor impedance: 0-5 k $\Omega$ ,
- Selectable language: GR, FR, ES, EN, DE, NL, IT, NR, PA, SE,
- Keypad on control panel in language: English,
- Included:
  - User manual in English(PDF),
  - Power cord,
  - ECG 3-lead cable,
  - Power supply,
- Dimensions: 11.4 x 6 x 3.1 cm,
- Weight: 0.5 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

