



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **WELCH ALLYN CP 200 Electrocardiograph ECG unit (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: through 4 new built-in 6.4 V / 6 Ah batteries,
- Frequency: 50/60 Hz,
- Welch Allyn CP 200 is a 12-lead resting ECG system that can also perform spirometry tests. The CP 200 is equipped with an alphanumeric keypad and LCD display. The display can show a quick color preview of the ECG reading along with the ability to edit settings. The Welch Allyn system provides full-page (8.25 x 11) ECG printouts and can store up to 50 ECG and 50 spirometry measurements in internal memory. The ECG system offers the ability to transfer patient data via SD card, UBS or your facility's EMR system,
- High-resolution thermal printer generates full-size standard and custom reports,
- Full alphanumeric keypad and LCD display for quick and easy entry of patient information,
- Advanced filters for optimal ECG recording quality ensure accurate readings,
- Configurable report formats allow information to be printed as desired,
- Battery powered for easy transport of the device to the patient,
- Instant-on function allows the device to be turned on quickly and perform ECG immediately,
- Technical specifications:
  - 12-lead resting ECG,
  - LCD screen for easy viewing,
  - Performs spirometry tests,
  - Stores up to 50 ECG readings and 50 spirometry readings,
  - Full-page printouts (8.25 x 11),
  - Configurable report formats,
- Welch Allyn cart:
  - Mobile on 4 casters (2 braking),
  - 2 side baskets for accessories,
  - Cable rack,
- Selectable menus in EN, DE, FR, ES, IT, PT, SE, NR, NL, FI, PL, RU, GR, HR, ML,
- Keypad on the control panel in English,
- Included:
  - User manual in English (PDF),
  - Power cord,
  - Mobile transport cart,
  - ECG 10-lead cable,
- Dimensions:
  - Welch Allyn CP 200: 41 x 40 x 16 cm,
  - Overall: 55 x 60 x 102 cm,
- Weight: 23 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

