



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



PHILIPS INTELLIVUE MP40 ANESTHESIA Vital function and anesthetic gas monitor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- PHILIPS MP40 ANESTHESIA patient monitors are portable monitors with measurement flexibility, making them meet the unique requirements for monitors when working in walk-in wards. Built on Philips' experience in patient monitoring, the IntelliVue monitor features flexible screen configuration, extensive clinical measurement menus, built-in transport facilities such as Event Surveillance, traditional 12-lead ECG recording with arrhythmia analysis and other useful capabilities. In keeping with Philips' adherence to providing customers with choice, IntelliVue monitors are compatible with a range of SpO2 measurement technologies and sensors manufactured by Philips, Masimo and Nellcor®. IntelliVue monitors are easy to use and can operate as part of a hospital-wide network platform. They can be configured according to the severity of the patient's condition, the standards of practice adopted in the department or the requirements of specific procedures,
- Technical specifications:
 - Panel screen, made with 12" TFT LCD technology and operating at SVGA resolution,
 - Convenient handle and sturdy case make it easy to carry,
 - A navigation dial makes it easy to enter information and navigate the screen,
 - IntelliVue technology - access information from various systems in the hospital network - pharmacy, laboratory, ECG management, radiology department PACS, server-based hospital systems and others - and deliver it to the patient monitor,
 - Audible and visual alarms,
 - Navigation joystick,
 - Ports on the monitor:
 - PS/2:
 - Mouse,
 - Keyboard,
 - Ethernet x2,
 - Lan x1,
- Philips M3001A module:
 - ECG /RESP - Brand new 3 - lead cable with extension,
 - SPO2 - Brand new silicone sensor,
 - NBP - Cuff with air tubing,
 - TEMP - Temperature connector socket, (no sensor included, can be retrofitted at additional cost),
- PHILIPS IntelliVue G5 M-1019A Gas Analysis Module:
 - Continuous real-time monitoring of respiratory gases and anesthetics (simultaneously),
 - IntelliVue G5 integrates completely with Philips IntelliVue MP40/50/60/70/80/90 anesthetic patient monitors,
 - Automatic identification of up to two anesthetic agents,
 - Inspiratory and expiratory levels of all gases,
 - Respiratory rate (awRR) calculated from the CO2 curve,
 - Minimum alveolar concentrations (MAC),
 - Long-term stability enables simple, infrequent calibration,
 - Individual alarm thresholds for each parameter,
 - Quick access to the most commonly used functions,
 - Selectable menu in language: ENG,
- Included:
 - User manual in English(PDF),
 - Power cord,
 - ECG /RESP - Brand new 3 - lead cable with extension,
 - SPO2 - Brand new silicone sensor,
 - NBP - Cuff with air tubing,
- Dimensions:
 - Monitor: 37 x 22 x 33 cm,
 - Gas module: 30 x 22 x 8.5 cm,
- Weight: 13 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

