



## **BERCHTOLD ELEKTROM 505 Surgical diathermy (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 40 W,
- Berchtold's EKTROTOM 505 is a compact and powerful HF surgical device for monopolar and bipolar surgical procedures,
- The ELEKTROTOM 505 high-frequency generator generates three different types of current:
  - Non-modulated high-frequency current for smooth cutting and contact coagulation,
  - Highly modulated high-frequency current for spray coagulation (fulguration) and cutting electrode coagulation,
  - Special low-resistance current suitable for bipolar coagulation,
- Technical specifications:
  - Nominal frequency: 500 kc/sec.
  - Modulation pulse frequency: 35 kc /s,
  - Observable number of cycles: 10/30 <sup>2</sup>,
- Included:
  - Foot manipulator,
  - ERBE neutral wire,
  - Bipolar cable,
  - ERBE monopolar cable - 2 buttons - 4 mm shaft,
  - Power cord,
- Dimensions: 25 x 23 x 8.5 cm,
- Weight: 3.8 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

