



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MELAG VACUKLAV 24 BL+ Autoclave (New)**

- New,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2100 W,
- Vacuklav 24 BL+ is a class B steam sterilizer and belongs to the group of devices with water connection. During its operation in a modern single water circulation system, it protects both the instruments and the autoclave,
- Vacuklav 24BL+ can be connected directly to a demineralizer, such as MELAdem 40 - it then works in automatic mode, or to an external tank that has been previously filled with water,
- Programs:
  - Quick program S: 134°C, sterilization time: 3.5 min / duration: 20 min / drying time: 5 min,
  - Quick program B: 134°C, sterilization time: 5.5 min. / duration: 32 min. / drying time: 10 min,
  - Universal program: 134°C, sterilization time: 5.5 min. / duration: 35 min. / drying time: 30 min,
  - Protective program: 121°C, sterilization time: 20.5 min. / duration: 48 min. / drying time: 20 min,
  - Prion program: 134°C, sterilization time: 20.5 min. / duration: 50 min. / drying time: 20 min,
- Technical specifications:
  - Fixed water and sewage connection,
  - Capacity: 29 L,
  - Carrying capacity: 7 KG of tools / 2.5 KG of textiles,
- Included:
  - Power cord,
  - 3 x tray,
  - 1x tray rack,
  - 2.5 m - water supply connection hose,
  - Drain hose,
  - Tray holder,
  - Attached side handles for transporting the autoclave,
- Dimensions:
  - Chambers: 60 x Ø25 cm,
  - Overall: 49.5 x 83.5 x 42.5 cm,
  - Pull-out trays: 54 x 19 cm,
- Weight: 53 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 12-months for domestic market (Poland),
  - 12-months for the international market,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

