



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ASPEL ASCARD A4 ECG unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Emergency power supply: via 1 new built-in 12 V 1.3 Ah battery,
- Frequency: 50/60 Hz,
- Power: 20VA,
- ASPEL ASCARD A4 unit for ECG signal analysis and interpretation,
- Specifications:
 - ECG signal: 12 standard leads,
 - Printout: 3, 6, 12 (after automatic analysis) - real-time channel (paper width: 112 mm),
 - Paper: heat-sensitive, dust-free with a width of 112 mm,
 - Recording speed: 5 / 10 / 25 / 50 mm/s,
 - Test sensitivity: 2.5 / 5 / 10 / 20 mm/mV,
 - Digital filters: 50 Hz and 35 Hz; Anti-dry filter,
 - Built-in LCD screen: 240 x 128 characters,
 - Foiled, touch and readable control panel,
 - Functional and alphanumeric keypad,
 - Signaling of bad connection of individual electrodes,
 - Output: RS-232,
- Selectable menu in the following languages: PL, EN, DE, FR, IT, ES, CR, TR, CZ, RU, SV, LT,
- Included:
 - User manual in Polish (PDF),
 - 10-lead ECG cable,
 - Power cord,
- Dimensions: 31 x 23 x 7 cm,
- Weight: 3 kg,
- Cat 02,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

