



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MAQUET ALPHAMAXX 1133.02 B3 Operating O/R table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: via 2 new built-in 12 V / 10 Ah batteries,
- Frequency: 50/60 Hz,
- MAQUET ALPHAMAXX is an operating table developed in close cooperation with physiotherapists, surgeons and nursing staff. It simplifies daily work in the operating room, supports patient safety and comfort. It can be easily modified according to the needs of any surgical discipline in just a few quick steps,
- New upholstery (memory visco) with medical approval in black,
- Adjustable in multiple planes by wireless remote control or through a panel on the column:
  - Height adjustment in the range: 60 - 105 cm,
  - Longitudinal movement in range: 22 cm,
  - Zeroing - return to the horizontal position of the tabletop,
  - Flex and Reflex - bending the tabletop in half in both directions,
  - Side tilt right and left in range: 20°,
  - Vertical tilt - Trendelenburg and Anti-Trendelenburg position: 30°,
  - Backrest section adjustment: -40 / +80°,
  - Electric brake,
- Manual adjustment:
  - Headrest: -45 / +90°,
  - Leg section - up/down -90°/10°,
  - Leg section - lateral retraction: 75°,
- Specifications:
  - 5 segments,
  - Electrohydraulic drive,
  - Universal table top,
  - X - RAY translucency,
  - Base housing made of durable plastic,
  - Equipped with rails for attaching additional accessories,
  - Table adapted for C-arm,
  - Mobile, on double casters,
- Includes:
  - English user manual (PDF),
  - Power cord,
  - Control remote control,
- Possibility to retrofit the table at an additional cost with other accessories from our offer, examples shown in the last photos,
- Dimensions:
  - Headrest: 50 x 30 x 6 cm,
  - Upper back support segment: 52 x 24 x 6 cm,
  - Middle seat-back segment: 52 x 94 x 6 cm,
  - Thigh rest segment: 24 x 12 x 6 cm,
  - Leg rest segment: 24 x 64 x 6 cm,
  - Overall length: 232 cm,
  - Overall width: 57 cm,
  - Base dimensions: 52 x 105 x 22 cm,
- Maximum load capacity: 450 kg (patient + accessories),
- Weight: 340 kg,
- Has an up-to-date inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

