



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KTMED KTPS-01 Pulse oximeter (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Korea, (2017),
- Power supply: 230 V,
- Emergency power supply: through 5 new built-in AA NiCd / NiMH 1.2 V batteries,
- Frequency: 50/60 Hz,
- KTMED KTPS-01 pulse oximeter provides stable patient monitoring which is applicable to neonates, adults and the elderly,
- The possibility of attaching the pulse oximeter to the rail through a bracket with a mount as well as the possibility of using it anywhere thanks to its small size,
- Advantages:
 - Easy to use,
 - Compact and portable design with light weight,
 - High accuracy and unreliability,
 - Audible alarms with adjustable limits,
 - Easy-to-read display,
- Included:
 - Power supply,
 - SPO2 cable,
- Dimensions: 11.5 x 5 x 16 cm,
- Weight: 0.5 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

