



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MediTop MT 1000 ENT unit (New)**

- New,
- Made in the Netherlands,
- MT1000 is a versatile treatment unit that can be used multifunctionally. The treatment unit can be used as a mobile or stationary model. It is ideal for smaller practices, clinics or acute care,
- The MT1000 comes standard with 5 parking handles equipped with suction cup, cold light, LED light source and/or compressed air. The cabinet has 2 instrument trays, 1 storage drawer and 1 disposal drawer for used instruments and waste,
- MT1000 is suitable for multidisciplinary spaces. From outpatient clinics to emergency rooms. Because the MT1000 is customizable, many options are possible. The device is mobile as standard. If the device is used in a fixed space, it can be equipped with an integrated microscope, but still need a versatile treatment unit with an integrated microscope,
- The device is also available in a stationary version, allowing integration of a microscope. This is especially useful if you have limited space
- The touchscreen control panel is intuitive to use and contributes to the user-friendly concept of the MT1000.
- Instruments can be set up ready for use and always at hand. Instruments placed in holders are automatically activated when removed and deactivated when reinserted,
- MT 1000 standard version:
  - 2 instrument trays and 5 holders
  - 2 instrument trays and 5 parking holders,
  - 1 large storage drawer,
  - Storage drawer for used instruments and waste,
- Possible integrated functions (at extra charge)
  - Suction system (various models),
  - Compressed air system (various models),
  - LED light source (various models),
  - Ear syringe with warm water,
  - Mirror heating (preheating and/or hot air),
  - Microscope with integrated camera (stationary version only),
  - HF surgery/coagulation,
- Optional accessories (at extra charge):
  - AQUATOP XS - Ear syringe with warm water,
  - LED MEDITOP - Headlamp,
  - MEDITOP - Comfort chair,
  - SPIGGLE & THEIS - instruments & disposables,
- Technical Data:
  - Connection voltage: 230V,
  - Maximum power: 2300VA,
- Dimensions: L. 1200 x D. 623 x H. 942 mm,
- Warranty: 12 months,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

