



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MediTop MT 750 ENT unit (New)

- New,
- Made in the Netherlands,
- The MT750 is a mobile and compact ENT unit. The MT750 offers various storage options for instruments and supplies for minor procedures. An ideal treatment module for use in multi-functional spaces,
- The MT750 comes standard with 3 holders for suction, compressed air and cold light, among others,
- With four casters, the MT750 can be easily moved. MT750 is sturdy and easy to use can be placed in any room,
- Details of MT750 unit construction (as shown in the picture below):
 - A - Unit compartment (devices can be stored in the unit compartment),
 - B - Handles (Instrument-related functions are automatically turned on and off by removing and reinserting the handle from the holder),
 - C - Instrument storage (Instrument trays are located in this compartment),
 - D - Waste storage,
 - E - Used instruments,
 - F - Electronics compartment (Behind the panel is the electronics compartment. The main switch is also located here),
 - G - Instrument storage compartment (Drawer G provides space for instruments and supplies),
 - H - Main switch (on the side of the unit) (The MT750 Basic can be turned on and off using the main switch),
 - J - Control panel,
- MT 750 standard version:
 - 2 instrument trays and 3 holders,
 - 1 large storage drawer,
 - Storage drawer for used instruments and waste,
- Possible integrated functions (at extra charge)
 - Suction system,
 - LED light source with cold light,
 - HF surgery/coagulation,
- Optional accessories (at extra charge):
 - AQUATOP XS - Warm water ear syringe,
 - LED MEDITOP - Headlamp,
 - MEDITOP - Comfort chair,
 - ZEISSA - Microscopes
 - SPIGGLE & THEIS - instruments & disposables,
- Technical Data:
 - Connection voltage: 230V,
 - Maximum power: 2300VA,
 - Dimensions: L. 991 x D. 575 x H. 935 mm,
- Warranty 12 months,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

