



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **AGA RE-PRO-LIFT RPL-1060/EE Proctology chair (New)**

- New,
- Made in Germany,
- The chair consists of two parts,
- It has electric height and pelvic tilt adjustment, backrest adjustment by gas spring,
- Ergonomically shaped seat and backrest cushions provide a very high level of comfort. Integrating a foam wedge in the back section eliminates the need for additional neck support. There is a cutout in the seat section,
- Standard equipment:
  - A pair of standard rails attached to the side of the seat,
  - Two holes (right/left side) for the colposcope holder,
  - Grounded socket on the back of the column,
  - Cutout in seat cushion,
  - Stainless steel tray with handle,
  - Paper roll holder attached to the back of the seat,
  - Foot switches,
- Technical specifications:
  - Frame: Lacquered,
  - Chair dimensions:
    - Length: 1150 mm,
    - Back/seat section: 800 / 350 mm,
    - Width: 580 mm (640 mm with standard rails),
  - Height adjustment: 650 - 850 mm (electric motor),
  - Upholstery segments: 2,
  - Thickness of upholstery in the seat section: 85 mm,
  - Backrest upholstery thickness: 85 - 135 mm,
  - Seat section adjustment: +30° (electric motor),
  - Backrest adjustment: 0°/+54° (gas spring),
  - Trendelenburg position: 30°
  - Maximum load capacity of 200 kg.
- Has a current Technical Passport issued,
- Warranty 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

