



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AGA RE-PRO-MED RP-1058/E Proctology chair (New)

- New,
- Made in Germany,
- The chair consists of two parts,
- It has electric adjustment of pelvic tilt, backrest adjustment by gas spring,
- Armchair with a fixed height of 865 mm,
- Ergonomically shaped seat and back cushions provide a very high level of comfort. The integration of a foam wedge in the back section eliminates the need for additional neck craning. There is a cutout in the seat section,
- The chair has been retrofitted with a folding footrest for easy chair entry,
- Standard equipment:
 - A pair of standard rails attached to the side of the seat,
 - Two holes (right/left side) for the colposcope holder,
 - Grounded socket on the back of the column,
 - Cutout in seat cushion,
 - Stainless steel tray with handle,
 - Paper roll holder attached to the back of the seat,
 - Foot switches,
 - Folding footrest,
 - Knee pad,
- Specifications:
 - Frame: Lacquered,
 - Seat dimensions:
 - Length: 1150 mm,
 - Back/seat section: 800 / 350 mm,
 - Width: 580 mm (640 mm with standard rails),
 - Height adjustment: fixed height 865 mm,
 - Upholstery segments: 2,
 - Upholstery thickness in the seat section: 85 mm,
 - Thickness of backrest upholstery: 85 - 135 mm,
 - Seat section adjustment: +30° (electric motor),
 - Backrest adjustment: 0°/+54° (gas spring),
 - Trendelenburg position: 30°
 - Maximum load: 200 kg.
 - Folding footrest: 490 x 320 x 250 mm (W x D x H),
- Has a current Technical Passport issued,
- 12-month warranty,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

