



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



HEBU HB8880 Oscillating plaster saw (Ex-Display)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany, (2022),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 250 W,
- HEBU HB8880 oscillating saw is extremely handy and easy to use for precise work. The electrically powered set allows the operator to quickly and safely remove synthetic and plaster dressings,
- The medical plaster cutting saw is a reliable tool that combines patient safety, precision cutting and convenient operation. Thanks to modern technological solutions, it provides invaluable support for medical professionals in the process of removing plaster dressings,
- The specially designed mechanism of the plaster cutting saw minimizes the risk of damaging the patient's skin during the procedure. The saw blades made of high-quality surgical steel guarantee precise plaster cuts. The speed adjustment mechanism allows the operator to adjust the cutting speed for different types of plaster,
- Technical specifications:
 - Speed adjustment: YES,
 - Motor speed: about 12000 - 21000 min,
 - Oscillation: 200-350 cycles/sec,
 - Noise level (Class A): 75 dB,
 - Vibration: Typically lower than 2.5 m/s²,
- Features:
 - Risk of patient injury during operation minimized,
 - Speed control,
 - Easy operation,
 - Easy disc replacement,
 - Removal of plaster and synthetic dressings,
 - Features soft-start function,
 - The motor is cooled by a cold air circulation system, which also protects against the deposition of gypsum dust and synthetic resin residue, extending its life,
- Included:
 - User manual in English, German, French, Spanish, Italian,
 - Built-in power cord,
 - 2 blades with diameters: 50 and 65 mm,
 - 2 pcs key 11 for fixing the blades,
- Dimensions of the device: 8 x 9 x 25 cm,
- Weight: 1.6 kg,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

