



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **SCHILLER CARDIOVIT AT-1 ECG unit (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Switzerland,
- Power supply: 230 V,
- Emergency power supply: via 1 new built-in 12 V / 1.3 Ah battery,
- Frequency: 50/60 Hz,
- Power: 28 VA,
- AT-1 allows operator to perform fast and accurate 12-lead ECG analysis in a three-channel format,
- With its compact form factor, very easy operation and built-in battery, it is a good device for all applications, which can include inpatient use in the office, at the bedside in the hospital and in outpatient settings,
- English language control panel keyboard,
- Rubberized buttons, easy to read control panel,
- Built-in thermal printer with a paper width of 90 mm (block or roll),
- Technical specifications:
  - The apparatus has two modes of operation:
    - Automatic mode,
    - Manual mode,
  - Recording speed: 5, 25, 50 mm/sec,
  - Test sensitivity: 5/10/20 mm/mV or automatic,
  - The instrument has three built-in filters:
    - Zero line filter: 0.05 (default), 0.15 or 0.30 Hz,
    - Network interference filter: 50, 60 Hz or off,
    - Muscle interference filter: 25 or 35 Hz on or off,
  - Test copy function,
- Included:
  - Power cord,
  - 10-lead ECG cable,
  - 4 pcs limb electrodes,
  - 6 pcs pre-cardiac electrodes,
- Device dimensions: 29 x 21 x 7 cm,
- Total weight: 3 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

