



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NEXUS Specialty bed for complex/trauma care (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the UK,
- Power supply: 230 V,
- Emergency power supply: through 4 built-in NaCa 6 V / 12 Ah batteries -replaceable at additional cost,
- Frequency: 50/60 Hz,
- Power: 1910 VA,
- Mattress (memory visco) with medical approval in blue,
- The NEXUS CB3SP is a complex care swivel bed designed specifically for the treatment and rehabilitation of patients after complex injuries, including patients with multi-organ injuries, spinal, neuro, pelvic and spinal cord injuries,
- The NEXUS CB3SP bed makes both the patient and clinician's use of complex critical care as safe and comfortable as possible, and is proven to reduce patient stress, recovery time and the risk of secondary complications,
- Technical specifications:
 - Electric mattress height adjustment: 59 - 95 cm,
 - Electric adjustment: footrest / knee 0° / 45°,
 - Electric adjustment of the side flaps of the mattress angle: 0° / 90°,
 - Electric adjustment of Trendelenburg / Anti-Trendelenburg position: 17°,
 - Electric adjustment of bilateral side tilt: 30°,
 - Electric backrest adjustment: 0° / 60°,
 - Electric adjustment via wired remote control:
 - Mattress height adjustment: 59 - 95 cm,
 - Trendelenburg / Anti-Trendelenburg position adjustment: 17°,
 - Adjustable side tilt: 30° on both sides,
 - Built-in CPR safety system,
 - Integrated digital weighing system,
 - Retractable side rails with 3-degree adjustable tilt along with the ability to retract them,
 - Built-in angle indicators,
 - SELF LEVEL button - highest / lowest position,
 - Built-in safety buttons,
 - Safety key system (full function / limited function),
 - Manual drip height adjustment within 72 cm,
- Mobile bed with central brake / guidance system,
- Included:
 - Manual in physical form,
 - Power cord,
 - Wired control remote control,
 - 2 pcs double-sided armrest,
 - 2 pcs armrest mounting bracket,
 - 2 pcs. security key:
 - Master - full function,
 - Default - limited function,
 - 3 segment footrest,
 - Drip rod mounting bracket,
 - 3 pcs traction support with guide roller,
 - Lower traction support guide wheel,
 - 2 pcs short traction support bracket,
 - 2 pcs long traction support bracket,
 - 2 pcs. head position stabilizer with overlay,
 - 2 pcs. head support bracket,
 - 6 pcs. rectangular bilateral side pads to brace the patient's body,
 - Metal operating strap,
 - Extendable drip tray with adjustable height,
 - Mobile storage compartment for accessories,
- Dimensions:
 - Upper mattress: 200 x 105 cm,
 - Lower mattress: 200 x 52 cm,
 - Overall: 220 x 120 x 70/105 cm,
- Maximum load capacity: 175 kg,
- Weight: 236 kg,
- It has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

