



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **SCHMITZ OPX MOBILIS RC 30 Operating O/R table (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: via 4 new built-in 5 V / 12 Ah batteries,
- Frequency: 50/60 Hz,
- Upholstery with medical approval in black (upholstery can be replaced with new upholstery at additional cost),
- SCHMITZ OPX RC 30 universal operating table with electro-hydraulic adjustment via remote control,
- All table positions can also be controlled manually with a hydraulic pump,
- Mobile table on 4 double covered wheels, 1 additional directional wheel, central brake, activated by foot switch,
- Technical specifications:
  - Electric adjustment by remote control:
    - Height adjustment: 72 - 108 cm,
    - Trendelenburg position: 30°,
    - Anti-Trendelenburg position: 30°,
    - Lateral tilt: -20 / 0 / 20°,
  - Manual adjustment:
    - Headrest tilt: -45 / 0 / 25°,
    - Backrest tilt: -50 / 0 / 70°,
    - Footrest inclination: -90 / 0 / 20°,
    - Footrests abduction angle: 55°,
  - Manual adjustment with a manual foot switch using a hydraulic pump:
    - Height adjustment: 72 - 108 cm,
    - Trendelenburg position: 30°,
    - Anti-Trendelenburg position: 30°,
    - Lateral tilt: -20 / 0 / 20°,
  - Table suitable for X-ray examination through the translucency of the table top,
- Included:
  - User manual in English (PDF),
  - Power cord,
  - Power supply,
  - Wired remote control,
- Upholstery has signs of wear that do not interfere with comfortable use,
- Dimensions:
  - Overall: 59 x 210 x 72/108 cm,
  - Headrest: 52 x 30 x 4 cm,
  - Backrest: 55 x 63 x 4.5 cm,
  - Seat: 52 x 48 x 4 cm,
  - Leg rest segment: 26.5 x 51 x 4 cm,
- Maximum load capacity:
  - Safe working load: 135 kg,
  - Safe working load in horizontal position: 185 kg,
- Weight: 210 kg,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

