



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PENTAX EG 2980K Videogastroscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- Technical specifications:
 - Viewing direction: forward view,
 - Maximum angle of field of view: 140°,
 - Depth of field: 4 - 100 mm,
 - Outer diameter: 9.8 mm,
 - Working channel diameter: 2.8 mm,
 - Working length: 1050 mm,
 - Total length: 1360 mm,
 - Max. Bend angle of the endoscope tip up and down: 210°/120°
 - Max. Right and left tip bending angle: 120°/120°,
- Included:
 - Original case,
 - Lid,
- Case dimensions: 62 x 42 x 14 cm,
- Weight: 4 kg,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

