



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **KLS MARTIN ML 301 (wall-mounted) Treatment and operating light (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 75 VA,
- Technical specifications:
  - Light intensity: 40,000 Lux,
  - Diameter of the illuminated field: 180 mm,
  - Canopy equipped with a halogen bulb 12 V / 50 W,
  - Lamp set on a two-part arm with a total length of 135 cm,
  - Ability to adjust the position of the canopy horizontally and vertically in the range of 360°,
  - Free height adjustment possible,
  - Built-in hand grip for free control of the lamp,
- Dimensions:
  - Canopy Ø: 34 cm,
  - Reflector diameter: 27.5 cm,
  - Total (lamp folded): 115 x 80 x 20 cm,
- Total weight: 8.5 kg,
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
  - 6-months for domestic market (Poland),
  - 3-months for the international market,
  - Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
  - The warranty period does not include bulbs and other light sources installed in the device,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

