



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



BAKMED. PF-222.4 Negatoscope (with stand) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 70 VA,
- Bulb: 11 W (2 pcs),
- X-ray image viewers, also called negatoscopes, are designed to evaluate X-ray images in operating rooms, doctors' offices, emergency rooms, X-ray laboratories, dental offices, etc. They are adapted to work as a stationary device (fixed to the wall) or portable (placed on a table),
- Screen luminance: 3500nt (+/- 20%),
- Mountable bottom stand to carry out a more stable working range,
- Included:
 - Power cord,
- Dimensions: 32 x 38 x 8 cm,
- Weight: 3.2 KG,
- Has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
 - The warranty period does not include bulbs and other light sources installed in the device,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

