



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EMALED 200 D PLUS LED treatment light (ceiling) (New)

- New,
- Made in Germany,
- A lamp for ceiling mounting,
- Advantages and features of the lamp:
 - Very suitable for examination and treatment rooms,
 - Light with high contrast and low energy consumption,
 - Shadeless illumination and good color reproduction,
 - Low energy consumption of test light of only 10-15 W,
 - Long LED life of at least 50,000 hours ensure low operating costs,
 - Low service costs and long service life ensure safe operation and low maintenance costs,
 - The light can be placed in any position thanks to the height-adjustable spring arm and 360 ° horizontal rotation,
 - High mobility and good positioning of the lamp body,
 - Natural color reproduction,
- Technical data:
 - Maximum light intensity (at 1 m): 70,000 lux,
 - Illumination intensity: 30 - 100 %,
 - Operating distance: 70 - 140 cm,
 - Light field diameter: 15 cm
 - Color temperature: 4 500 K,
 - Color reproduction index 95 Ra,
 - Led lifetime 50,000 h,
 - Canopy diameter 24cm,
 - 7 LEDs,
 - Fixed focal length,
 - Main voltage: 230 V,
 - Frequency 50 - 60 Hz,
 - Power of the lamp head 10VA,
 - Weight of the lamp head 1 kg,
 - Total weight 12 kg,
- Power supply and ceiling mounting kit included,
- Suspension tube of 680mm (it is possible to order this model with suspension tube of 1500mm at an additional cost),
- Has a current Technical Passport issued,
- 24-month warranty,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

