



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **EMALED 200 M PLUS LED treatment lamp (mobile with battery) (New)**

- New,
- Made in Germany,
- Mobile lamp (mobile base with battery),
- Advantages and features of the lamp:
  - The EMALED 200M PLUS mobile lamp features an integrated emergency power supply and can therefore be used for up to 24 hours without a power connection,
  - Very suitable for examination and treatment rooms,
  - High mobility and good body positioning,
  - Light with high contrast and low energy consumption,
  - No moving parts in the lamp head,
  - Shadeless illumination and good color reproduction,
  - Low energy consumption of test light of only 10-15 W,
  - Long LED life of at least 50,000 hours ensure low operating costs,
  - Low service costs and long service life ensure safe operation and low maintenance costs,
  - Powder coated aluminum housing,
  - High mobility and good lamp body positioning,
  - Barely noticeable heating of the lamp head,
  - Natural color reproduction,
  - 3-stage brightness adjustment,
  - Two autoclavable holders are included (additional holders can be ordered at any time),
- Specifications:
  - Maximum light intensity (at 1m): 70,000 lux,
  - Illumination intensity: 30 - 100 %,
  - Operating distance: 70 - 140 cm,
  - Light field diameter: 14 cm,
  - Color temperature: 4 500 K,
  - Color reproduction index 95 Ra,
  - Led lifetime 50,000 h,
  - Canopy diameter 24cm,
  - 7 LEDs,
  - Fixed focal length,
  - Main voltage: 230 V,
  - Frequency 50 - 60 Hz,
  - Head power
  - Power of the lamp head 10VA,
  - Weight of the lamp head 1 kg,
  - Total weight 37 kg,
- Includes mobile base with 230V power supply, Swivel base, spring arm, charger with battery and EMALED 200 lamp head,
- Has a current Technical Passport issued,
- Warranty 24 months,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

