



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MPW-53 Laboratory centrifuge (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 60 W,
- MPW 53 is a small laboratory centrifuge, its design ensures ease of use, safe operation and wide application in medicine, biochemistry and other laboratories where there is a need to separate a small amount of the test preparation into components of different densities under the influence of centrifugal force,
- Advantages:
 - Small size and simple control system,
 - Choice of operating time up to 15 min or continuous operation,
 - Choice of two centrifugation speeds,
 - Optimization of acceleration and deceleration time,
- Technical specifications:
 - Maximum acceleration: 3120 x g,
 - 2 RPM modes:
 - 3500 RPM,
 - 5800 RPM,
 - Spinning time: 0 - 15 min,
- Included:
 - Power cord,
 - Swivel rotor 11101 - standard rotor for centrifugation of blood tubes or 6 x 15/10 ml conical tubes (tube outer ø: 18 mm),
- Dimensions: 22 x 27 x 18 cm,
- Weight: 4 KG,
- Has up-to-date inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

