



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FRESENIUS KABI INJECTOMAT Syringe infusion pump (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in France,
- Power supply: 230 V,
- Emergency power supply: through 5 new built-in 6 V / 1.8 Ah batteries,
- Frequency: 50/60 Hz,
- Power: 15 W,
- FRESENIUS KABI INJECTOMAT electronically controlled single syringe pump for cooperation with central power supply and data management system,
- Pump calibrated to work with syringes of different types and volumes and different manufacturers,
- Features:
 - Versatile syringe pump with dose calculation,
 - For administering drug treatment,
 - Advanced functions, including dose calculation in units of weight, programmable bolus, initiating dose,
 - Easy to use,
 - Intuitive operation,
 - Clear, interactive menus,
 - Robust and proven design,
 - Library of 50 medications,
 - Visual and audible alarms,
- Technical specifications:
 - Range, dosing rate in range: from 0.1 to 200 ml/h (increments of 0.1 ml/h) (up to 1200 ml/h, programmable),
 - Infusion volume: from 0.1 to 999.9 ml (increments of 0.1 ml),
 - Manual bolus: 50 to 1200 ml/h. (increments of 50 ml/h),
 - Blue monochrome graphic LCD display, 70 mm x 35 mm,
 - Automatic syringe recognition: 5 / 10 / 20 / 30-35 / 50 / 60 ml,
- The pump supports syringes:
 - BD PERFUSOR 50 ml,
 - BD PLASTIPAK 50/30/20/10/5 ml,
 - BD PLASTIPAK WWD 50/20/20 ml,
 - BRAUN OMNIFIX 50 ml,
 - BRAUN PERFUSOR 50/20 ml,
 - CODAN 50/10 ml,
 - DIDACTIC LINE 50 ml,
 - DIDACTIC PERFUSION 50 ml,
 - DISPOMED INFUSECT 50 ml,
 - DISPOMED PERFUSECT 50 ml,
 - FRESENIUS INJECTOMAT 50/10 ml,
 - FRESENIUS P-SPRITLE 50 ml,
 - ICOGAMMA PLUS 50 ml,
 - SANPOL 50 ml,
 - KABIFILL 50 ml,
 - MAROOMED 50 ml,
 - MOWOJECT 50/30/20/10/5 ml,
 - PENTA 60/30/20/10 ml,
 - PIC INDOLOR 50 ml,
 - POLFA 50/20/5 ml,
 - TERUMO 50/30/20/10/5 ml,
 - TUTOJECT TYPET 50 ml,
 - BRAUN OMNIFIX 30/20 ml,
 - FRESENIUS MED.CARE 30 ml,
- Selectable menus in PL, FR, EN, DE, NR, SE, DU, IT, ES, GR, RU, HU, CS, SK, TR, RU, PT, BR,
- Included:
 - User manual in English (PDF),
 - Power cord,
- Dimensions: 34.5 x 13.5 x 16 cm,
- Weight: 2.1 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for a period of 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!
In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

