



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS CF EVIS EXERA II H180AL Video colonoscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Japan,
- Technical specifications:
 - Channel width: 3.7 mm,
 - Working length: 168 cm,
 - Field of view: 170,
 - Depth of field: 2-100,
 - Outer diameter: 12.8 mm,
 - Maximum up/down angle: 180 degrees,
 - Maximum right/left angle: 160 degrees,
- Included:
 - Suction valve,
 - Water/air valve,
 - Stopper (plug) for endoscope cleaning,
 - Original carrying case,
 - User manual in English (PDF),
- It has a current inspection and is ready for operation,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

