



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Double-sided X-ray protective gown with collar (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- X-ray radiation protection gown is designed to reconcile comfort and uncompromising secondary radiation protection,
- Made of lightweight multilayer and flexible lead rubber for excellent performance,
- Protects the front, back, shoulders and thyroid of the wearer during x-rays,
- Fastens with Velcro,
- Dimensions:
 - Chest width: 53 cm,
 - Hip width: 57 cm,
 - Overall length: 120 cm,
- Weight: 8 KG,
- Has a current inspection and is ready for use,
- Issued Technical Passport (Service Report) valid for 12 months,
- Warranty:
 - 6-months for domestic market (Poland),
 - 3-months for the international market,
 - Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
- Financing options (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

